




Created by  
King's College London  
Biomedical Engineering  
Students



# Medical Devices Made Simple



**DON'T MISS THE  
KIDS' CORNER**  
A special fun page  
just for kids

# Glucometer

A device that quickly measures blood sugar (glucose) levels

## Key Features

**Quick and Easy:**  
Provides results in just a few seconds

**Portable and lightweight:** small enough to fit in a pocket or bag

**User-Friendly:**  
Simple Interface, making it easy for anyone to use



## Physical Characteristics

**Consists of a:**  
**Glucometer** – measures glucose levels

**Disposable test strips** – to insert the blood into the device

**Lancet** – a small needle to prick the finger to retrieve a drop of blood

### Interesting Facts

- 1 Before glucometers, doctors tested urine instead of blood to estimate sugar levels, but it wasn't as accurate or immediate.
- 2 The first glucometer, called the Ames Reflectance Meter, was invented in the 1970s and was originally the size of a large book!
- 3 Dogs can be trained to detect low blood sugar in people with diabetes by sensing changes in their owner's scent—sometimes even before a glucometer picks it up!

### Why It's Important

Knowing your blood sugar (glucose) levels is essential for managing diabetes and maintaining overall health.

### Where is it used?

Ideal for hospitals, clinics, and home monitoring - helps to detect high and low blood sugar levels - essential for people with diabetes to manage to condition

### What does it measure

Blood sugar levels is measured in milligrams per decilitre (mg/dL)

## How to Use the Device

**Step 1:** Insert the test strip into the glucometer

**Step 2:** Use the lancet to prick your fingertip for a drop of blood

**Step 3:** Place the blood drop on the test strip.

**Step 4:** Wait a few seconds for the device to display your glucose level

# Hearing Aids

A compact, non-invasive device that help individuals with hearing loss regain access to sounds

## Key Features

Provides **clear sound amplification** instantly, enhancing hearing with minimal effort

**Easily worn** throughout the day and discreetly carried in your pocket or bag

Designed to be **simple and comfortable** for users of all ages, with easy-to-use controls



## Physical Characteristics

**Size:** Compact and discreet, designed to fit comfortably in or behind the ear

**Weight:** Lightweight (around 2 – 6 grams)

**Power:** Battery-operated (with long-lasting life, many featuring rechargeable options)

### Where is it used?

Ideal for individuals with hearing loss/disorder, to improve communication and sound clarity.

### Interesting Facts

① The first hearing aid was invented in the late 19<sup>th</sup> century, often in the form of a large, cumbersome device.

② Modern hearing aids can be customized with different features, including Bluetooth connectivity for phone calls or streaming music.

③ Many hearing aids are now digital, offering noise reduction, feedback cancellation, and directional microphones for better clarity

### What is it capable of

1. **Speech clarity:** Improves understanding of speech, especially in noisy environments.
2. **Environmental Sounds:** Enhances awareness of surrounding sound, such as traffic or household noises.

### Why It's Important

Hearing aids help individuals with hearing loss regain access to sounds, improving communication, social interaction, and overall QOL. They can reduce the risks associated with hearing loss, such as isolation, cognitive decline and depression.

## How to Use the Device

**Step 1:** Place the hearing aid in/behind your ear

**Step 2:** Turn the device on using the power button

**Step 3:** Adjust the volume as needed

**Step 4:** Turn off to save battery life

# Inhaler

An inhaler is a handheld device that delivers medicine as a spray for easier breathing

## Key Features

**Quick and Easy** –  
Delivers medicine in seconds

**Portable and Lightweight** – Inhalers fit inside your pocket to keep with you all day

**Effective Relief** –  
Directly targets the lungs for fast action



## Physical Characteristics

**Size:** pocket-sized

**Weight:** Lightweight (around 100 grams)

**Different colours indicate their purpose:**

**Blue:** Reliever Inhaler

**Orange/Brown:** Prevent Inhaler

**Green/Purple/Red:**

Combination Inhaler

### How does it work?

Inhalers deliver medicine directly to the lungs, relaxing airways and reducing swelling

Some inhalers use pressurised gas, while others are powder based

### Interesting Facts

- 1 The first modern inhaler was invented in 1956 by George Maison for his daughter
- 2 Some inhalers use spacers to make it easier to inhale the medicine
- 3 Cats can get asthma and can be treated using an inhaler similar to that of humans.

### Types of Inhalers:

- **Reliever Inhalers:** used to relieve asthma symptoms when they occur
- **Prevent Inhalers:** contain steroids to prevent asthma symptoms from occurring – used everyday
- **Combination:** combine the effects of both reliever and preventer inhalers – taken everyday

### Why It's Important

An inhaler helps people with asthma, allergies, and lung conditions breathe more easily by delivering medicine straight to the lungs

## How to Use the Device

**Step 1:** Shake the inhaler before use

**Step 2:** Exhale fully, then place in your mouth

**Step 3:** Press down on the inhaler while breathing in

**Step 4:** Hold your breath for a few seconds, then breathe out slowly

# Nebuliser


A device that turns liquid medicine into a fine mist that can be inhaled through a mouthpiece

## Key Features

**Effective Delivery** – Ensures medication reaches deep into the lungs for optimal treatment

**User-Friendly** – Suitable for patients who struggle with inhalers

**Versatile** – can be used with a variety of medication!



## Interesting Facts

- 1 Jean Sale-Girons invented the first nebulizer, the 'pulverisateur' in 1858 – it used a handpump to create mist
- 2 NASA has researched nebulisers to ensure astronauts can receive medication in space
- 3 Not just for Humans! Nebulisers are also used for animals, including horses, cats, dogs and birds

## Why it's important

A nebulizer is helpful for a variety of conditions:

**Chronic Obstructive Pulmonary Disease (COPD)**

**Severe Asthma Attacks**

**Bronchiectasis and Pulmonary**

**fibrosis** – use with saltwater solution to manage phlegm build up

## Types of Nebulisers:

**Jet Nebulisers** – use compressed air to turn medication into mist (most common)

**Ultrasonic Nebulisers** – use high frequency vibrations (mainly for hospital use)

## How to Use the Device

**Step 1:** Fit all the components together and fill the medicine cup

**Step 2:** Sit upright to ensure effective delivery of the medication

**Step 3:** Turn on the compressor and breathe in the mist slowly (10–15 minutes)

**Step 4:** After each use, clean it as instructed to keep it hygienic and working well



# Portable BP

A compact, non-invasive device that quickly measures blood pressure (BP) and Pulse Rate

## Key Features

### Quick and Easy:

Provides accurate blood pressure readings in just a few seconds

### Portable and lightweight:

Easily carried in your pocket or bag for convenient monitoring on the go

### User-Friendly:

Simple interface designed for anyone to use without prior training



## Physical Characteristics

**Size:** Approximately 4"x3"

**Weight:** Lightweight (around 200 – 300 grams)

**Power:** Battery-operated (with auto power-off to save battery)

### Where is it used?

Ideal for hospitals, clinics, and home monitoring— useful for people managing hypertension needing regular blood pressure tracking.

### Interesting Facts

- ① The first blood pressure monitor was invented in 1881 by Samuel Siegfried Karl Ritter von Basch
- ② High blood pressure (hypertension) is often called the "silent killer" because it has no symptoms
- ③ Some smart blood pressure monitors can sync with mobile apps for easy tracking and analysis

### What does it measure

1. **Systolic Pressure** (top number): The pressure in your arteries when your heart beats.
2. **Diastolic Pressure** (bottom number): The pressure in your arteries between beats.
3. **Pulse Rate:** The number of heart beats per minute.

### Why It's Important

Monitoring your blood pressure regularly helps in detecting hypertension early, preventing heart disease, stroke, and other cardiovascular issues.

## How to Use the Device

**Step 1:** Wrap the cuff around your upper arm

**Step 2:** Press the power button to start

**Step 3:** Remain still while the cuff in/deflates

**Step 4:** Wait and view the result

# Pulse Oximeter

A compact, non-invasive device that quickly measures blood oxygen saturation (SpO<sub>2</sub>) and pulse rate

## Key Features

**Quick and Easy:**  
Provides accurate readings in just a few seconds

**Portable and lightweight:** Easily carried in your pocket or bag

**User-Friendly:**  
Simple Interface designed for anyone to use without prior training



## Physical Characteristics

**Size:** Approximately 1"x3"

**Weight:** Lightweight (around 100 – 150 grams)

**Power:** Battery-operated (with auto power-off to save battery)

### Where is it used?

Ideal for hospitals, clinics, and home monitoring—especially useful for anyone wanting a quick check on their vital signs.

### Interesting Facts

- ① Pulse oximeters were originally developed for aviation to help pilots monitor oxygen levels at high altitudes
- ② They have become essential tools for early detection of respiratory issues, including monitoring COVID-19 complications
- ③ Their accuracy can be influenced by factors such as nail polish or poor circulation in cold fingers

### What does it measure

A pulse oximeter uses light to detect the difference in light absorption between oxygenated and deoxygenated blood. It then calculates the **oxygen saturation (SpO<sub>2</sub>)** and **pulse rate** based on the changes in light sensed in 5 seconds

### Why It's Important

A pulse Oximeter helps monitor your blood oxygen levels and heart rate which are essential for your overall respiratory and cardiovascular health

## How to Use the Device

**Step 1:** Turn on the pulse oximeter

**Step 2:** Place your fingertip into the sensor tip

**Step 3:** Wait a few seconds until the reading

**Step 4:** Displays your oxygen level and pulse

# Digital Thermometer

A medical device used to measure body temperature accurately

## Key Features

### Fast and Accurate:

Provides temperature readings in seconds with high precision

### Multiple

#### Measurement modes:

can be used orally, rectally, or underarm, offering flexibility

### Safe and Hygienic:

Comes with disposable probe covers to prevent cross-contamination.

## Infrared Thermometer

Non-contact devices that measure temperature from the forehead or ear using infrared technology

– Ideal for quick hygienic, and large-scale screenings



## Handheld Thermometer

Compact, Portable, and typically used for oral, axillary, or rectal measurements. These are suitable for home and clinical use

### Normal Temperature Ranges:

- ✓ Oral: 36.1°C – 37.5°C
- ✓ Rectal: 36.6°C – 38.0°C

### Temperature Display Options:

- ✓ **Celsius (°C):** Standard international unit.
- ✓ **Fahrenheit (°F):** Commonly used in the United States.
- ✓ **Kelvin (K):** More commonly used in medical settings, primarily scientific applications.

### Interesting Facts

- 1 In 1593, Galileo Galilei invented a simple water thermometer, allowing people to measure temperature changes for the first time.
- 2 Mercury thermometers used a shiny liquid that moved up and down with temperature changes, but now we often use digital ones because they're safer and easier to read.
- 3 Animals like snakes and some bugs can sense temperature changes without a thermometer, using special heat-sensing organs!

## How to Use the Device

**Step 1:** Turn on the thermometer

**Step 2:** Place the probe in your mouth or under your armpit

**Step 3:** Wait a few seconds until the reading stabilises

**Step 4:** Displays your temperature



# Learn More About Medical Devices

Click on the links to go to websites that can give you more information

## GLUCOMETER

<https://www.nhs.uk/conditions/cgm-and-hcl-for-diabetes/>

<https://www.diabetes.org.uk/about-diabetes/looking-after-diabetes/technology/flash-glucose-monitors-and-continuous-glucose-monitors>

## BLOOD PRESSURE MONITOR

<https://www.bhf.org.uk/information-support/heart-matters-magazine/medical/tests/blood-pressure-measuring-at-home>

<https://www.england.nhs.uk/ourwork/clinical-policy/cvd/home-blood-pressure-monitoring/>

## HEARING AIDS

<https://www.thcp.co.uk/articles/>

<https://www.nhs.uk/conditions/hearing-aids-and-implants/>

<https://www.nidcd.nih.gov/health/hearing-aids>

## NEBULISER

<https://www.asthmaandlung.org.uk/symptoms-tests-treatments/treatments/nebulisers>

<https://www.lhch.nhs.uk/nebulisers>

<https://www.pari.com/uk/nebulisers/>

## INHALER

<https://www.asthmaandlung.org.uk/healthcare-professionals/adult-asthma/choosing-inhaler-device>

<https://www.rightbreathe.com/>

<https://www.plymouthhospitals.nhs.uk/display-pil/>

## PULSE OXIMETER

<https://www.england.nhs.uk/wp-content/uploads/sites/52/2022/02/pulse-oximeter-easy-read-2022-digital.pdf>

<https://www.cuh.nhs.uk/patient-information/how-to-use-a-pulse-oximeter/>

## THERMOMETER

<https://thermometer.co.uk/content/295-thermometry-basics>

<https://education.nationalgeographic.org/resource/thermometer/>

<https://www.sciencing.com/infrared-thermometers-work-4965130/>